



USAID
FROM THE AMERICAN PEOPLE

Chest X-ray Taking Procedures Training for X-ray Technicians/ Radiographer

“Conventional Imaging”

Ms. KHIN YADANAR MOE

Consultant (TB CXR Training), IDDS Project/ Myanmar

Content



X-ray Film



Film Construction



Types of Film Based on Their Application



Film Artifacts



Film Handling and Storage



Intensifying Screen



Factors Affecting Image Quality



Film Cassette



Film Processing



Limitations

Film – Screen Radiography

- In film-screen radiography, a sheet of film with a light-sensitive emulsion on both sides is **sandwiched** between two intensifying screens
- In order to overcome the thicker screen, sandwiching the double-sided film between two screens serves to add the optical densities of the two emulsions together

Reference: John Ball , Tony Price, Chesney's Radiographic

X-ray Film

Reference: John Ball , Tony Price, Chesney's Radiographic

X-ray Film

- A media that makes a permanent record of the image
- To record the image on the film
- Images are stored as a latent image

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Construction

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Construction

I) Polyester plastic base:

- Support layer for fragile emulsion
- Must be clear, strong and consistent thickness (0.18 mm)
- Tinted pale blue or blue-gray (reduces eye strain)
- Coated on one or two sides with emulsion
 - *single emulsion film (better detail)
 - *double emulsion film (less detail)

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Construction – Cont.

2) Adhesive Layer

- To provide uniform surface over which the emulsion can be coated uniformly

3) Substratum layer

- To keep emulsion layer and base adhering to each other during the coating stage and processing

Reference: John Ball , Tony Price, Chesney's Radiographic


Film Construction – Cont.

4) Emulsion layer:

- Acts as photographic active layer
 - *activated by light and radiation to create image
- Consists of a mixture of gelatin solution and silver halide crystal
 - *typical emulsion consists of silver bromide (98%) and silver iodide (2%)
 - *crystal size – 1.0 to 1.5 μ in diameter

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Construction – Cont.



Film Construction

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Construction – Cont.

“Latent Image Formation”

- Exposure of silver – iodo – bromide grains to light photons emitted by the screen
- Affects of grain size and distribution:
 - The bigger the average grain size, the higher the **speed** of the film
 - The more grain distribution, the lower the **contrast**
 - The bigger the crystal, the higher the **graininess** (clumping of the crystal)
 - Latent image is made visible by chemical processing

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Construction – Cont.

“White and Dark Appearance on Film”

- **Dark** – silver halide crystals exposed to photons
*after processing, it turns to black metallic silver
- **White** – no crystals exposed
*silver halide is wash away by processing

Reference: John Ball , Tony Price, Chesney’s Radiographic

Film Construction – Cont.

5) Supercoat

- To protect the layer of **gelatin**
- To reduce damage from scratches and pressure

Gelatin

- Used as a suspending medium and binding agent for the silver halide particles
- Can be easily spread on the film base when warm
- Firmly exists on the base as a gel when cool
- Flexible and does not crack easily on bending

Reference: John Ball , Tony Price, Chesney's Radiographic

Types of Film based on their Application

Reference: John Ball , Tony Price, Chesney's Radiographic

Types of Film based on their application

1) Screen Films:

- Most commonly used
- Sensitive to blue light emitted by intensifying screens and direct actions of x-rays
- Used in cassettes with intensifying screens
- High speed

Reference: John Ball , Tony Price, Chesney's Radiographic

Types of Film based on their application – Cont.

2) Non-screen / Direct exposure films:

- Has thicker coat of emulsion used without an intensifying screen
- Depends mainly on action of x-ray
- Four times faster than that of screen films
- Must be manually processed because of thick emulsion
- Uses :
 - a) To detect intra-ocular foreign bodies
 - b) In dental with intra-oral cardboard

Types of Film based on their application – Cont.

3) Mammography film:

- Single coated
- Designed to be used with single intensifying screen
- Combination must be fast to deliver minimum dose to the glandular tissue

Reference: John Ball , Tony Price, Chesney's Radiographic

Types of Film based on their application – Cont.

4) Duplicating Film:

- Used to copy radiograph
- Original cassette to be copied is inserted into a cassette whose opaque front has been replaced by a pane of clear glass
- Special duplicating film is placed with emulsion side down onto the radiograph

Reference: John Ball , Tony Price, Chesney's Radiographic

Types of Film based on their application – Cont.

Differences between single coated and double coated x-ray film		
Characteristic	Single side	Double coated
Radiation dose	Increased	Decreased
Detail	Increased	Decreased
Parallex effect	No	Yes
Contrast	Decreased	Increased

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Artifacts

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Technical artifacts (caused by X-rays tube)

No	Artifacts	Causes	Remedy
1	Over exposure	<ul style="list-style-type: none">- High exposure factors- High exposure on fast speed screens	<ul style="list-style-type: none">- Use optimum factors- Use film – screen combination
2	Under exposure	<ul style="list-style-type: none">- Low exposure factors- Low exposure on slow speed screens	<ul style="list-style-type: none">- Use optimum factors- Use film – screen combination

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Technical artifacts (caused by X-rays tube) – Cont.

No	Artifacts	Causes	Remedy
3	Double exposure	- Single film is exposed twice	- Always separate exposed and unexposed cassettes
4	Film fog / noise	- Use of damaged cracked cassette	- Carefully use the cassette - Replace the damaged cassette on time

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Technical artifacts (caused by X-rays tube) – Cont.

No	Artifacts	Causes	Remedy
5	Grid cut off effect	<ul style="list-style-type: none">- Due to inaccurate position of patient on the table- Grid is not completely inserted into the table	<ul style="list-style-type: none">- The patient's midsagittal plane should be aligned with the mid axis of the table- Grid should be completely inserted into the table

Reference: John Ball , Tony Price, Chesney's Radiographic

2. Positioning artifacts (caused by position of patient, cassette and tube)

No	Artifacts	Causes	Remedy
1	Image blurring	- Due to improper respiration	- Ask patient to practice inspiration, expiration and holding of breath
2	Marker effect	- Due to the use of wrong marker or side placement	- Use accurate markers and on the proper side

Reference: John Ball , Tony Price, Chesney's Radiographic

2. Positioning artifacts (caused by position of patient, cassette and tube) – Cont.

No	Artifacts	Causes	Remedy
3	Movement unsharpness	- Due to patient movement during exposure	<ul style="list-style-type: none">- Explain the procedure to the patient- Use immobilization devices- Reduce exposure
4	Metallic Artifact	- Presence of metallic object between tube and film	<ul style="list-style-type: none">- Ask the patient to remove the radiopaque objects

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Processing artifacts

No	Artifacts	Causes	Remedy
1	Film fog	<ul style="list-style-type: none">- Damaged safe light- Accidental exposure of white light or any general light	<ul style="list-style-type: none">- Check safe light periodically or within 15 days- Cassette should be closed- Door interlocked
2	Age fog	<ul style="list-style-type: none">- Use of expired film	<ul style="list-style-type: none">- Always use film before the expiry date and check the date regularly

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Processing artifacts – Cont.

No	Artifacts	Causes	Remedy
	Age fog	- Storing of film in high temperature and high humidity conditions	- Air conditioner should be installed and check properly
3	Cresenteric effect	- Due to image of finger on film during handling	- Always handle the film carefully with the help of two fingers during loading and unloading

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Processing artifacts – Cont.

No	Artifacts	Causes	Remedy
4	Stain effect (various colours may appear on film during processing)		
	Brown colour stain	<ul style="list-style-type: none">- Improper fixing- Exhausted fixer- Improper washing	<ul style="list-style-type: none">- Proper fixing- Change fixer on time- Proper washing

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Processing artifacts – Cont.

Artifacts	Causes	Remedy
Yellow colour stain	<ul style="list-style-type: none">- Oxidized developer- Use of old developer	<ul style="list-style-type: none">- Developer tank should be covered and preservative added- Change developer solution on time
Static elastic effect (stripe)	<ul style="list-style-type: none">- Due to friction between two films	<ul style="list-style-type: none">- Film should be handled properly and film box should be placed vertically

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Handling & Storage

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Handling

1. Do not flex
2. Hands must be clean
3. Film is sensitive to pressure, scratches, light, x-rays, heat, moisture, electricity and age

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Storage

1. Proper film storage is required for **high quality and long-lasting images**
2. The location must be clean, dry and light tight
3. 40-60% humidity and 69° F or 10 – 24 ° C temperature
4. Avoid storing near chemical fumes which can fog the film
5. Safe from radiation exposure

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Storage

6. Expiration date clearly visible
7. Should be stored on edge (like books in a library)
8. Do not stack the film boxes horizontally because film in the bottom may show pressure artifacts
9. Always used older film first

Reference: John Ball , Tony Price, Chesney's Radiographic

Intensifying Screen

Reference: John Ball , Tony Price, Chesney's Radiographic

Intensifying Screen

- An image amplifier **converting the aerial image** which is relatively made of a few x-rays photons into many thousand times of light photons
- Converts incident **x-ray photons to visible light** which then exposes the silver halide emulsion on the film
- The image is formed **by processing the latent image** using dark room procedures

Reference: John Ball , Tony Price, Chesney's Radiographic

Intensifying Screen Construction

- Polyester plastic base – support layer
- Phosphor layer – active layer (X-rays photons converted to light photons) *Photoelectric effect
- Reflective layer – increases screen efficiency
- Protective coating

Reference: John Ball , Tony Price, Chesney's Radiographic

Intensifying Screen Construction – Cont.

Intensifying Screen Phosphor;

1. Rare Earth phosphor (emits green light) **Standard**
 - most efficient and commonly used
2. Calcium Tungstate (emits blue light)
 - not efficient

Reference: John Ball , Tony Price, Chesney's Radiographic

Intensifying Screen Construction – Cont.

I. Rare Earth Screen:

- Higher DQE (detective quantum efficiency)
**the percentage of x-rays absorbed by the screen*
- Higher x-rays absorption abilities
- Higher CE (Conversion efficiency)
**the amount of light emitted for each x-rays absorbed*
- More light emitted per x-rays absorbed by the screen

Reference: John Ball , Tony Price, Chesney's Radiographic

Types of Intensifying Screen with different speed

Screen

High resolution

Regular or Standard

Fast

Speed

Slow

Medium

Fast

Reference: John Ball , Tony Price, Chesney's Radiographic

Types of Intensifying Screen with different speed

- Screen speed is defined as a relative number to demonstrate how efficiently x-rays are converted into useable light
 - *Screen speed ranges from 100 (slow) to 1200 (fast)*
 - *Routine - 200 to 800, high detail – 50 to 100*

Faster screen speed - reduces patient exposure

Quantum mottle

Image detail ↓

Image noise (speckled background on the image) ↑

Reference: John Ball , Tony Price, Chesney's Radiographic

Advantages and Disadvantages of Intensifying Screen

Advantages	Disadvantages
A screen can absorb 20 – 40 times more x-rays than film alone	Less detail than direct exposure
Reduced exposure	Formation of quantum mottle and unsharpness
Reduced patient dose	Reduced detail
Increased x – ray tube life	Maintenance of screen

Reference: John Ball , Tony Price, Chesney's Radiographic

Advantages and Disadvantages of Intensifying Screen – Cont.

Advantages	Disadvantages
Contrast good	Screen artifact
Prevent motion lack of sharpness	-
Suitable for pediatric, geriatric and large body parts	-
Short developing time	-

Reference: John Ball , Tony Price, Chesney's Radiographic

Factors Affecting Image Quality

Reference: John Ball , Tony Price, Chesney's Radiographic

Factors affecting image quality

- Rare Earth screen should be used to be a good image quality.
- The use of intensifying screens lowers spatial resolution (good differentiation between two nearby objects) compared with direct-exposure radiographs

Reference: John Ball , Tony Price, Chesney's Radiographic

Factors affecting image quality – Cont.

- Spatial Resolution:
 - The higher the lp/mm the smaller the object that can be imaged
 - Very fast screens - 7 lp/mm
 - Fine-detail screens - 10 lp/mm
 - Direct-exposure screens - 50 lp/mm

*Spatial resolution is expressed by the number of **line pairs per millimeter** (lp/mm)*

Reference: John Ball , Tony Price, Chesney's Radiographic

Factors affecting image quality – Cont.

- Noise reduces the image contrast
 - *caused by high kvp and fast screens used*
- Artifacts (unwanted information in the image) reduces the image quality
 - Small scratches and dirty screen

*Screens should be **cleaned once each month** with manufacturer's cleaner with antistatic compounds*

Reference: John Ball , Tony Price, Chesney's Radiographic

Film Cassette

Reference: John Ball , Tony Price, Chesney's Radiographic

What is a Film Cassette?

“A film cassette is a container for exposed or unexposed film.”

Functions:

- To hold intensifying screen and protect them from damage
- To exclude all light from entering the cassette and fogging the film
- To maintain a close and uniform contact between the film and screens
- To exclude dust and dirt from the sensitive screens

Reference: John Ball , Tony Price, Chesney's Radiographic

Care of cassette

- Do not scratch the surface
- Keep clean (use mild soap and water on a regular basis)
- Do not allow the screen to become wet
- Record the date of cleaning

Reference: John Ball , Tony Price, Chesney's Radiographic

Care of cassette – Cont.

- Hold it gently and store in a standing position
- Do not store cassette near sources of heat
- Do not leave the cassette open
- Ensure the screen is fully dry before reloading the cassette

Reference: John Ball , Tony Price, Chesney's Radiographic

Cassette storage

- Upright and away from radiation
- Do not stack
- Should always be loaded and ready to use
- Avoid humidity and dust

**Number the cassettes so repeat problems can be easily identified*

Reference: John Ball , Tony Price, Chesney's Radiographic



Cassette

Consultant's own training material

Film Processing

Reference: John Ball , Tony Price, Chesney's Radiographic

Manual Film Processing:

- Involves the processing of the film by chemicals with the help of a person
- Latent image is converted into visible image in this process
- This process contains (5) steps:
 - 1) Developing
 - 2) Rinsing
 - 3) Fixing
 - 4) Washing
 - 5) Drying

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Developing

Developing

Is a chemical process in which the latent image is converted into a visible image

To convert metallic silver into black metallic silver by reduction process

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Developing – Cont.

- ***A developer solution contains:***
 - i. Solvent
 - ii. Developing agent
 - iii. Activator
 - iv. Preservative
 - v. Restrainer

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Developing – Cont.

i. Solvent

- Commonly used solvent is **water**
- The solvent **dissolves chemicals** and also aids **ionization**

ii. Developing agent

- Hydroquinone, phenidone and metol
- The purpose of developing agent is **to convert exposed AgBr crystals to black metallic silver**
- Hydroquinone is responsible for high contrast
- Metol responsible for grey shades

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Developing – Cont.

- Hydroquinone is always used in combination with metol to **shorten the developing time**
- Works more efficiently when the solution temperature is less than 20 degree Celsius

iii. Activator

- Commonly used hydroxide, sodium carbonate and sodium metaborate
- Used to open the pores of the film and allow the developing agent to do their work
- Optimum pH is range from 10 to 11

Reference: John Ball , Tony Price, Chesney's Radiographic

I. Developing – Cont.

iv. Preservative

- Use sodium sulphate or potassium metabisulphite
- **Decrease the rate** of oxidation of hydroquinone (**developing agent**)
- Increases the **life of the developer solution**

v. Restrainer

- Use potassium bromide solution
- It is also called the anti fogging agent

Reference: John Ball , Tony Price, Chesney's Radiographic

2. Rinsing

- When the x-ray film is removed from the developer, some chemicals remain on the film
- The film should be rinsed to remove these chemicals
- To stop the reaction of development of the developer and neutralization the basicity of the residual developer solution
- In a water bath, rinsing the film for **30 seconds**

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Fixing

- Is the process of **removing the unexposed AgBr** without damaging the image formed by metallic silver
- It also hardens the gelatin emulsion
- The optimum temperature is 18-24 degree Celsius
- Fixing time is **one to four minutes**

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Fixing – Cont.

- ***The fixing solution contains:***
 - i. Solvent
 - ii. Acidifier
 - iii. Clearing agents
 - iv. Hardening agent
 - v. Preservative

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Fixing – Cont.

i. Solvent

- Water is normally used as solvent
- Used to dissolve the chemicals

ii. Acidifier

- Commonly use sulphuric or acetic acid
- Neutralize the basicity of the residual developer which remains on the film
- It also provides the suitable medium for the fixer and hardener to act

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Fixing – Cont.

iii. Clearing agents

- Commonly use sodium or ammonium thiosulphate salt
- Also acts as fixing agent

iv. Hardening agent

- Commonly used chemical is potassium Alum
- To hard the gelatin emulsion
- Decrease the physical injury of the film

Reference: John Ball , Tony Price, Chesney's Radiographic

3. Fixing – Cont.

v. *Preservative*

- Commonly use chemical sodium sulphite
- Increase the life of the fixer solution
- **Protects** the fixing agent from **decomposition** (damage)

Reference: John Ball , Tony Price, Chesney's Radiographic

4. Washing

Washing:

- After fixing the film, the film must be washed with water
- Washing **removes the residual** processing solution and fixing chemicals
- If these chemical are **not removed** the image will **discolor** and fade
- Normally, x-ray film should be washed with **distilled water or tap water**
- Washing time is normally about **20 min** at 20 degree Celsius

Reference: John Ball , Tony Price, Chesney's Radiographic

Limitations

Reference: John Ball , Tony Price, Chesney's Radiographic

Limitations in Conventional Radiography

- The radiographic speed is fixed and not possible to adjust patient dose
- Narrow exposure latitude (low visualization of soft tissue and bone)
- Fixed brightness and grey-scale that cannot be adjusted
- Many toxic chemicals are used

Reference: John Ball , Tony Price, Chesney's Radiographic

Limitations in Conventional Radiography – Cont.

- High repeat exposure rate
- Imaging archiving is difficult
- Time intensive
- Increased radiation dose

Reference: John Ball , Tony Price, Chesney's Radiographic

THANK YOU !